


LAPIPTM

PROTOCOL

The LAPIP[®] Protocol for Saving Ailing & Failing Implants



The LAPIP protocol is a minimally-invasive approach to treating ailing and failing dental implants.

The LAPIP protocol utilizes the unique PerioLase[®] MVP-7TM, a Nd:YAG free running pulsed laser requiring specific knowledge and skills derivative of the LANAP protocol.

To you, the patient, this means:

- Protect implant investment
- No cutting
- No bone grafting
- No stitches
- Less pain
- Less tissue removed
- Less down time
- Less recovery time
- Less bleeding
- Less swelling

The LAPIP protocol:

- Allows removal of diseased tissue without removing any of the healthy tissue
- Eliminates the infection
- Helps maintain the height of the tissue around implants
- Minimizes pain and discomfort to the patient
- Closes the periodontal pocket wound to allow healing to take place

The LAPIP protocol is a no-risk first line treatment that preserves other treatment options. Controlling peri-implant disease makes it possible to save and protect your implant investment.

Treatment Procedure

- A) Periodontal probe indicates excessive pocket depth
- B) PerioLase[®] MVP-7TM laser energy selectively vaporizes bacteria, endotoxins, and diseased tissues; and removes titanium corrosion contaminants in soft tissue around implant
- C) Ultrasonic scaler tips remove debris from implant
- D) Bone is modified to release stem cells
- E) PerioLase MVP-7 finishes cleaning pocket, stimulates natural growth factors, and aids in sealing pocket closed so new germs cannot enter
- F) Healing of gums to clean surface occurs
- G) Bite trauma is adjusted
- H) Regeneration of bone around implant

First Last DDS

Street Address

City, ST zip

(800) 555-1212

www.yourwebsite.com

